

Bones Muscles Joints

The Spectra Knee CPM

kinetec
USA

An unbeatable combination of good looks and mechanical innovation sets Spectra apart as the Knee CPM of choice. The entire spectrum of needs for knee patients can now be addressed with a carriage that accommodates all patients, age 8 to 80 – without any additional accessories. Lightweight and exceptionally easy to carry, Spectra offers a quick set-up and practical operational modes that enhance patient comfort, compliance and results.

Features:

- Long-life DC motor and screw-drive mechanism for years of reliable performance
- Anatomically correct alignment increases patient compliance and decreases joint stress
- Wide carriage (10 inches) accommodates the smallest children, biggest athletes; offers a complete range of treatment
- INTERNALLY contained medical grade transformer eliminates handling extra components
- Exclusive Load Reverse ensures patient protection when excess force is exerted on the joint
- Lightweight, only 24 pounds!
- Sleek design incorporates carry handles – from either side
- Includes comfortable, easy-to-apply patient pad kit
- Patient lock-out prevents accidental program changes and potential injury
- Waterproof digital hand control snaps easily into foot support
- Pause setting allows 15 minutes of rest at flexion or extension
- 2-year full warranty ensures reliability
- UL/C-UL assures compliance to the highest standards



Technical Specifications

Range of Motion:	Hyperextension (-10 °) to full Flexion (120 °)
Weight:	24 pounds (53 kg)
Length:	37 inches (94 cm)
Full leg size:	28 to 41 inches (71 to 104 cm)
Tibia:	15 to 23 inches (38 to 58 cm)
Femur:	13 to 18 inches (33 to 46 cm)

NEW

Spectra Essential Knee CPM.

"Easy-to-use" and "intuitive" describes the Spectra Essential. The new remote control allows for easy adjustments of the R.O.M., Speed and Pauses in extension and in flexion. The large, clear display gets an immediate feedback which is always very useful for the motivation of the patient. All other features are the same as the Spectra Knee CPM except the treatment modes are not available on the Essential.

0815-27-514

